



Reichert SPR Small Molecule Analysis

95 Daltons and 201 Daltons Analyte Surface Plasmon Resonance (SPR) Analyses Using the Reichert SR7500DC and Carboxymethyl Dextran Slides

This application note showcases both the very low noise of the SR7500DC Dual Channel SPR system and the utility of Reichert's carboxymethyl dextran slide, CMD500k, for use with small molecule interactions. The interaction between two small molecule inhibitors, methanesulfonamide (95 daltons) and 4-carboxybenzenesulfonamide (4-CBS) (201.2 daltons), and the enzyme carbonic anhydrase (CA II) is studied. Since carboxymethyl dextran sensor chips were first introduced into the marketplace about 20 years ago, they have become the industry standard for obtaining SPR data. They are used to reduce non-specific binding and increase sensitivity. The capture layer of the Reichert CMD500k sensor chip (catalog number 13206066) is a 500,000 dalton carboxymethyl dextran hydrogel coupled to a self-assembled monolayer formed on the gold surface. The combination of the SR7500DC System and the CMD 500k sensor chip yields excellent results for these interactions.

Experimental

The experimental conditions are summarized in the following table:

Ligand	Analyte	Analyte Concentrations	Association Time	Dissociation Time
CA II	Methanesulfonamide	1066, 533, 266, 133 and 66.6 μM	1 min	1 min
CA II	4-CBS	5, 2.5, 1.25, 0.625, 0.3125 and 0.156 μM	1 min	1.5 min

Results




FIGURE 1

Target 1: Methanesulfonamide (95 Da) binding to CA II

About 6,000 mRIU of carbonic anhydrase II was amine coupled to the dextran surface and then the binding of each inhibitor was followed over a series of concentrations (see table). The binding of the low molecular weight inhibitor, methanesulfonamide, generates small but very detectable responses owing to the high sensitivity of the SR7500DC. The kinetic fit from Scrubber (Biologic Software) yields an equilibrium dissociation constant (K_D) of 652 μ M, while the Langmuir isotherm (inset) yields a K_D of 410 μ M. Both values are within the range of values previously published by other researchers.




FIGURE 2

Target 2: 4-CBS (201 Da) binding to CA II

4-carboxybenzenesulfonamide binding to carbonic anhydrase II for a series of concentrations (see table) is shown here. Each concentration was injected at least in duplicate. The values obtained for the equilibrium dissociation constant (K_D) are 809 nM for the kinetic fit from Scrubber (Biologic Software) and 940 nM for the Langmuir isotherm (inset). Both are within the range of values published in the literature for this binding pair using SPR and isothermal titration calorimetry (ITC).

Learn more about how Reichert pushes the limits of detection and sensitivity in label-free interaction. To improve the quality of your protein interaction research, visit www.ReichertSPR.com

®Trademark of Reichert Technologies Company Notice: No freedom from infringement of any patent owned by Reichert or others is to be inferred. Because use conditions and applicable laws may differ from one location to another and may change with time, Customer is responsible for determining whether products and the information in this document are appropriate for Customer's use and for ensuring that Customer's workplace and disposal practices are in compliance with applicable laws and other government enactments. The product shown in this literature may not be available for sale and/or available in all geographies where Reichert is represented. The claims made may not have been approved for use in all countries. Reichert assumes no obligation or liability for the information in this document. References to "Reichert" or the "Company" mean the Reichert legal entity selling the products to Customer unless otherwise expressly noted. NO WARRANTIES ARE GIVEN; ALL IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ARE EXPRESSLY EXCLUDED.

Corporate Office
 Reichert Technologies
 3362 Walden Avenue
 Buffalo, New York 14043, USA
 Tel: +1 716-686-4500
 Fax: +1 716-686-4555
 Toll Free USA: 1-888-849-8955
 reichertspr.lifesciences@ametec.com

European Service Center
 Carl-von-Linde Str. 42-85716
 Unterschleissheim/Munich, Germany
 Tel: +49 89 315 8911 0
 Fax: +49 89 315 891

www.reichertspr.com